

Rabbit Anti-GGGGS Linker (R4E06) mAb

Catalog No.	RGP16001A RGP16001B	Quantity:	100 µg 1.0 mg
Alternate Names:	G4S, G4S link, G4S Linker, GGGGS, GS Linker		
Description:	Affinity purified rabbit monoclonal (clone R4E06) that targets the GGGGS Linker.		
Specificity:	GGGGS Linker		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Rabbit IgG		
Formulation:	Sterile-filtered PBS containing 0.09% sodium azide		
Clone:	R4E06		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FC: 0.25 µg per 10 ⁶ cells in 100 µl ELISA: suggested dilution 1:5,000-1:10,000 IF: suggested dilution 1:2,000-1:8,000 IHC: suggested dilution 1:40,000-1:160,000 Western blot: suggested dilution 1:1,000-1:6,000		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

