

Rabbit Anti-SDMA/Symmetric Dimethylarginine (SAb2291) mAb

Catalog No.	RGP16301A RGP16301B	Quantity:	100 µg 1.0 mg
Alternate Names:	Symmetric Dimethylarginine, SDMA		
Description:	Affinity purified rabbit monoclonal (clone SAb2291) that targets Symmetric Dimethylarginine for use in ELISA		
CAS Number:	30344-16-4		
Specificity:	Symmetric Dimethylarginine (SDMA)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Rabbit IgG		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	SAb2291		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

