

GLRA3

Anti-Glycine receptor Alpha3 (N424/83) mAb

Catalog No.	RHB14103A RHB14103B	Quantity:	100 µg 1.0 mg
Alternate Names:	Glycine receptor subunit alpha-3, GLRA3		
Description:	Monoclonal Mouse Anti-Glycine receptor Alpha3 (N424/83) is purified by affinity chromatography from cell culture supernatant.		
UniProt ID (target):	O75311		
Target:	Glycine receptor subunit alpha-3, GLRA3		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	N424/83		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	IHC, Western blot		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

