

## IL2RA

### Mouse Anti-Human CD25/IL2RA (BA9) mAb

<b>Catalog No.</b>	RHB95801A RHB95801B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, TAC antigen, p55, CD25		
<b>Description:</b>	Monoclonal Mouse Anti-Human CD25/IL2RA (BA9) targets Human Interleukin-2 receptor subunit alpha and is purified by affinity chromatography.		
<b>UniProt ID (target):</b>	P01589		
<b>Target:</b>	Human Interleukin-2 receptor subunit alpha		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	Sterile filtered 0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	BA9		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Applications:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.