

## APP

### Mouse Anti-Human APP/Amyloid beta (10D5)

<b>Catalog No.</b>	RHC12504A RHC12504B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Amyloid-beta precursor protein, ABPP, APP, Alzheimer disease amyloid A4 protein homolog, Alzheimer disease amyloid protein, Amyloid precursor protein		
<b>Description:</b>	Mouse monoclonal antibody (clone 10D5) which targets Human Amyloid beta.		
<b>UniProt ID (target):</b>	P05067		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Human Amyloid beta		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Clone:</b>	10D5		
<b>Purification:</b>	Protein A or G affinity chromatography		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	ELISA, SPR		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.