

CEACAM5

Research Grade Mitumomab

Catalog No.	RHC21001A RHC21001B	Quantity:	100 µg 1.0 mg
Alternate Names:	CEA-Scan, IMMU-4, CAS: 154361-48-5		
Description:	Research Grade Mitumomab is a murine monoclonal antibody which targets Human Carcinoembryonic antigen-related cell adhesion molecule 5 (CEA, CD66e)		
UniProt ID (target):	P06731		
Specificity:	Human Carcinoembryonic antigen-related cell adhesion molecule 5 (CEA, CD66e)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgG2a		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Mitumomab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

