

PTHLH

Mouse Anti-Human PTHrP/PTHLH (1-34) (SAA0326) mAb

Catalog No.	RHC93101A RHC93101B	Quantity:	100 µg 1.0 mg
Alternate Names:	Parathyroid hormone-related protein, PTH-rP, PTHrP, Parathyroid hormone-like protein, PLP		
Description:	Affinity purified Mouse anti-human PTHrP/PTHLH (1-34) clone SAA0326 is a monoclonal antibody that targets both human and mouse Parathyroid hormone-related protein (1-34)		
UniProt ID (target):	P12272		
Specificity:	Human, Mouse Parathyroid hormone-related protein (1-34)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG1		
Formulation:	Sterile filtered 0.01M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	SAA0326		
Cross Reactivity:	Mouse Parathyroid hormone-related protein (1-34)		
Applications:	Neutralization, Western blot		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

