

## MMP8

### Rabbit Anti-Human/Mouse/Rat MMP8 (R1S18) mAb

<b>Catalog No.</b>	RHD55102A RHD55102B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Neutrophil collagenase, Matrix metalloproteinase-8, MMP-8, PMNL collagenase, PMNL-CL		
<b>Description:</b>	Monoclonal rabbit anti-human/mouse/rat MMP8 (R1S18) was affinity purified from cell culture supernatant.		
<b>UniProt ID (target):</b>	P22894		
<b>Target:</b>	Human, Mouse, Rat MMP8		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Rabbit IgG		
<b>Formulation:</b>	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium azide		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A/G purified from cell culture supernatant		
<b>Clone:</b>	R1S18		
<b>Applications:</b>	<b>FCM:</b> suggested dilution 1:20-1:100 <b>IF:</b> suggested dilution 1:50-1:200 <b>IHC:</b> suggested dilution 1:100-1:200 <b>IP:</b> suggested dilution 1:20-1:50 <b>Western blot:</b> suggested dilution 1:1000-1:2000		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)