

GLP1R

Anti-Human GLP1R (SAA0519) mAb

Catalog No.	RHE40101A RHE40101B	Quantity:	100 µg 1.0 mg
Alternate Names:	Glucagon-like peptide 1 receptor, GLP-1 receptor, GLP-1-R, GLP-1R		
Description:	Mouse monoclonal Anti-Human GLP1R (SAA0519) was affinity purified from cell culture supernatant.		
UniProt ID (target):	P43220		
Specificity:	Human Glucagon-like peptide 1 receptor		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2a kappa		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	SAA0519		
Purification:	Protein A/G purified from cell culture supernatant		
Applications:	ELISA, IF, IHC		
Cross Reactivity:	Monkey		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

