

MMP14

Anti-Human MMP14 (LEM-2/15) mAb

Catalog No.	RHE75501A RHE75501B	Quantity:	100 µg 1.0 mg
Alternate Names:	Matrix metalloproteinase-14, MMP-14, Membrane-type matrix metalloproteinase 1, MT-MMP 1, Membrane-type-1 matrix metalloproteinase, MT1-MMP, MMP-X1		
Description:	Affinity purified mouse monoclonal antibody that targets human Matrix metalloproteinase-14		
UniProt ID (target):	P50281		
Specificity:	Human Matrix metalloproteinase-14		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG1 kappa		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	LEM-2/15		
Purification:	Protein A/G purified from cell culture supernatant		
Applications:	ELISA, ICC, IF, FC, Western blot EC ₅₀ = 78.15, clone LEM-2/15 binds with human MMP-14		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

