

SMAD2

Mouse Anti-Human Phospho-Smad2 (pT8) (SAA0289) mAb

Catalog No.	RHH20002A RHH20002B	Quantity:	100 µg 1.0 mg
Alternate Names:	Mothers against decapentaplegic homolog 2, SMAD 2, SMAD2, JV18-1, Mad-related protein 2, hMAD-2, MADR2, Mothers against DPP homolog 2, SMAD family member 2, MADH2, Smad2, hSMAD2, MAD homolog 2		
Description:	Mouse monoclonal antibody that targets human Mothers against decapentaplegic homolog 2 (SMAD 2)		
UniProt ID (target):	Q15796		
Specificity:	Human MAD homolog 2 (SMAD 2)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2a kappa		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	IHC, WB, IP		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

