

LYPD1

Mouse Anti-Human LYPD1 (SAA1492) mAb

Catalog No.	RHJ21001A RHJ21001B	Quantity:	100 µg 1.0 mg
Alternate Names:	Ly6/PLAUR domain-containing protein 1, LYPD1, Putative HeLa tumor suppressor, PHTS		
Description:	Mouse Anti-Human LYPD1 Antibody (SAA1492) is specific for human Ly6/PLAUR domain-containing protein 1		
UniProt ID (target):	Q8N2G4		
Target:	Human Ly6/PLAUR domain-containing protein 1 (LYPD1)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Clone:	SAA1492		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com