


## Angptl4

### Anti-Mouse ANGPTL4 (6G11) mAb

<b>Catalog No.</b>	RHJ51403A RHJ51403B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Angiopietin-related protein 4, Angiopietin-like protein 4 Hepatic fibrinogen/angiopietin-related protein, HFARP, Fasting-induced adipose factor, Secreted protein Bk89		
<b>Description:</b>	Neutralizing mouse monoclonal antibody (clone 6G11) which targets the C-terminus of Mouse Angiopietin-related protein 4.		
<b>UniProt ID (target):</b>	Q9Z1P8		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Mouse Angiopietin-related protein 4 (C-terminus)		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Clone:</b>	6G11		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Purity:</b>	> 95% by SDS-PAGE and SEC-HPLC		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Endotoxin:</b>	< 0.01 EU/µg by LAL analysis		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Manufacturer:</b>	Antibody System		
<b>Applications:</b>	Functional studies, neutralization, ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

#### SEC-HPLC



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)