


Mouse Anti-Semaglutide (GLP-1 analog) (SAb2275) mAb

Catalog No.	RHK13902A RHK13902B	Quantity:	100 µg 1.0 mg
Alternate Names:	LY3298176, GIP/GLP-1 RA, TZP		
Description:	Rabbit Anti-Tirzepatide polyclonal antibody purified by antigen affinity chromatography. Tirzepatide (LY3298176, GIP/GLP-1 RA, TZP) is a dual GIP/GLP-1 receptor agonist.		
UniProt ID (target):	P01275		
Specificity:	Human GLP-1, Semaglutide, Tirzepatide		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	SAb2275		
Applications:	ELISA		
Manufacturer:	Antibody System		
Storage & Stability:	Store at 2-8°C for immediate use, within 2 weeks. Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

Anti-Semaglutide Antibody (SAb2275) binds with Semaglutide, GLP-1, and Tirzepatide



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

SEC-HPLC

