

## TTR

### Mouse Anti-Human Transthyretin amyloidosis (6C1) mAb

<b>Catalog No.</b>	RHK25503A RHK25503B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Transthyretin, Prealbumin, TBPA, TTR, PALB, ATTR		
<b>Description:</b>	Mouse Anti-human transthyretin amyloidosis clone 6C1 is affinity purified neutralizing monoclonal antibody.		
<b>UniProt ID (target):</b>	P02766		
<b>Specificity:</b>	Human Transthyretin		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	Sterile filtered 0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A or G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	6C1		
<b>Applications:</b>	ELISA, Neutralization, Western blot, SPR		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**