

scFv hlgG1 Fc(C220A/L234A/L235A) Isotype Control (HyHEL-10) mAb

Catalog No.	RHK32401A RHK32401B	Quantity:	100 µg 1.0 mg
Description:	Affinity purified Human scFv hlgG1 Fc(C220A/L234A/L235A) recombinant monoclonal antibody (HyHEL-10) is used as an isotype control in flow cytometry, and <i>in vivo</i> control.		
Source:	Mammalian cell line		
Specificity:	Human IgG1 Fc(C220A/L234A/L235A) Fragment		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	HyHEL-10		
Extinction Coefficient:	1.63		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Molecular Weight:	103 kDa		
Concentration:	1.0 mg/ml		
Manufacturer:	Antibody System		
Applications:	FCM, Isotype Control, <i>In vivo</i> Control		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.