

ROBO1

Anti-Human ROBO1 (B5209B) mAb

Catalog No.	RHK45002A RHK45002B	Quantity:	100 µg 1.0 mg
Alternate Names:	Roundabout homolog 1, ROBO1, Deleted in U twenty twenty, H-Robo-1		
Description:	Human monoclonal antibody that targets human Roundabout homolog 1 was expressed in a mammalian cell line and purified using Protein A/G affinity chromatography.		
UniProt ID (target):	Q9Y6N7		
Target:	Human Roundabout homolog 1		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	B5209B		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.