

## TRBV9

### Mouse Anti-Human TRBV9 (SAA1760) mAb

<b>Catalog No.</b>	RHK86201A RHK86201B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	T cell receptor beta variable 9, TRBV9		
<b>Description:</b>	Murine monoclonal antibody that targets human T cell receptor beta variable 9 was expressed in a mammalian cell line and purified using Protein A/G affinity chromatography.		
<b>UniProt ID (target):</b>	A0A0B4J1U6		
<b>Target:</b>	Human T cell receptor beta variable 9		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cell line		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	SAA1760		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

