

HLA-A, Tax

Mouse Anti-Human HLA-A2-peptide (LLFGYPVYV) Complex (T3F2) mAb

Catalog No.	RHM03018A RHM03018B	Quantity:	100 µg 1.0 mg
Alternate Names:	HLA class I histocompatibility antigen, A alpha chain, HLA-A, Human leukocyte antigen A, HLAA, Protein Tax-1, Protein X-LOR, Protein PX, Trans-activating transcriptional regulatory protein of HTLV-1		
Description:	Mouse Anti-Human HLA-A2-peptide (LLFGYPVYV) Complex (T3F2) monoclonal antibody is affinity purified from cell culture supernatant.		
UniProt ID (target):	P04439 & P03409		
Target:	Human HLA-A2-peptide LLFGYPVYV (Protein Tax-1)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG1		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	T3F2		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com