

## SEZ6L2

### Anti-Human SEZ6L2 (SAA1851) mAb

<b>Catalog No.</b>	RHN01601A RHN01601B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Seizure 6-like protein 2, SEZ6L2, PSK		
<b>Description:</b>	Mouse monoclonal anti-Human Seizure 6-like protein 2 (SAA1851) was affinity purified from cell culture supernatant.		
<b>UniProt ID (target):</b>	Q6UXD5		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Human Seizure 6-like protein 2,		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Clone:</b>	SAA1851		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	PBS containing 50% glycerol, 0.5% BSA and 0.09% sodium azide		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for frequent use or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

