

SEC61A2

Anti-Human SEC61A2 (R2S74) mAb

Catalog No.	RHN86701A RHN86701B	Quantity:	50 µg 100 µg
Alternate Names:	Protein transport protein Sec61 subunit alpha isoform 2, Sec61 alpha-2		
Description:	Rabbit monoclonal anti-human Protein transport protein Sec61 subunit alpha isoform 2 (SEC61A2) is affinity purified from cell culture supernatant.		
UniProt ID (target):	Q9H9S3		
Source:	Mammalian cells		
Specificity:	Human Protein transport protein Sec61 subunit alpha isoform 2 (SEC61A2) Cross reacts with mouse and rat		
Isotype:	Rabbit IgG		
Clone:	R2S74		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Purity:	> 95% by SDS-PAGE		
Formulation:	TBS containing 0.05% BSA, 50% glycerol, 0.05% sodium azide		
Concentration:	1.0 mg/ml, lot specific		
Applications:	IF: 1:50 - 1:200 IHC: 1:50 - 1:100 IP: 1:20 Western blot: 1:500 - 1:1000		
Storage & Stability:	Store at 2-8°C for frequent use or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com