

Ifng

Rat Anti-Murine IFNG/IFN-gamma (XMG1.2) mAb

Catalog No.	RMB95401A RMB95401B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interferon gamma, IFN-gamma		
Description:	Recombinant monoclonal antibody rat anti-mouse IFNG/IFN-gamma (XMG1.2) was expressed in a mammalian cell line and purified using Protein A/G affinity chromatography.		
UniProt ID (target):	P01580		
Target:	Murine IFN-gamma		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Rat IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	XMG1.2		
Purification:	Protein A/G affinity chromatography		
Applications:	Flow cytometry		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com