

ApoE

Mouse Anti-Murine APOE (HJ6.3) mAb

Catalog No.	RMB98601A RMB98601B	Quantity:	100 µg 1.0 mg
Alternate Names:	Apolipoprotein E, Apo-E		
Description:	Recombinant antibody mouse anti-murine Apolipoprotein E (HJ6.3) was expressed in a mammalian cell line and purified using Protein A/G affinity chromatography.		
UniProt ID (target):	P08226		
Target:	Mouse Apo-E		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	Block, IHC		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.