

Fgf2

Anti-Mouse FGF2/bFGF (GAL-F2) mAb

Catalog No.	RMC38201A RMC38201B	Quantity:	100 µg 1.0 mg
Alternate Names:	Fibroblast growth factor 2, FGF-2, Basic fibroblast growth factor, bFGF, Heparin-binding growth factor 2, HBGF-2		
Description:	Affinity purified mouse monoclonal antibody that targets murine fibroblast growth factor 2. Useful for blocking and flow cytometry.		
UniProt ID (target):	P15655		
Specificity:	Mouse Fibroblast growth factor 2		
Cross Reactivity:	Human Fibroblast growth factor 2		
Isotype:	Mouse IgG1, kappa		
Clone:	GAL-F2		
Formulation:	0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Purity:	>95% as determined by SDS-PAGE		
Applications:	Blocking, FCM		
Storage & Stability:	Store at 2-8°C for immediate use, or prepare aliquots and store at -20°C to -80°C for up to 1 year. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

