

Mouse IgG Format Isotype Control for Bispecific Antibody (SAb2315)

Catalog No.	RMK81501A RMK81501B	Quantity:	100 µg 1.0 mg
Description:	Mouse IgG Format Isotype Control for Bispecific Antibody (SAb2315) used as an isotype control and <i>in vivo</i> control.		
Source:	Mammalian cell line		
Isotype:	Mouse IgG like (KiH)		
Clone:	Sab2315		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Applications:	In vivo Control, Isotype Control, Anti-Mouse CD3 & Anti-Lysozyme (HyHEL-10) Isotype Control		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

