

LMP1

Human Anti-EBV/HHV-4 LMP1/BNLF1 (B2A4) mAb

Catalog No.	RVV07005A RVV07005B	Quantity:	100 µg 1.0 mg
Alternate Names:	Latent membrane protein 1, LMP-1, Protein p63, BNLF-1		
Description:	Human monoclonal (clone B2A4) to Latent membrane protein 1 from Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)		
UniProt ID (target):	P03230		
Specificity:	Latent membrane protein 1 from Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)		
Source:	Mammalian cell line		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	B2A4		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for immediate use. Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com