

G

Mouse Anti-HeV/NiV Glycoprotein G (HENV-26) mAb

| | | | |
|---------------------------------|---|------------------|------------------|
| Catalog No. | RVV08003A RVV08003B | Quantity: | 100 µg 1.0 mg |
| Alternate Names: | Glycoprotein G, Attachment glycoprotein, Receptor binding protein, RBP, Nipah Virus (NiV) / Hendra Virus (HeV) | | |
| Description: | Murine neutralizing monoclonal to Hendra virus (isolate Horse/Australia/Hendra/1994) & Nipah virus Glycoprotein G (clone HENV-26) | | |
| UniProt ID (target): | O89343 & Q9IH62 | | |
| Specificity: | HeV & NiV Glycoprotein G | | |
| Source: | Mammalian cell line | | |
| Isotype: | Mouse IgG1 kappa | | |
| Formulation: | Sterile-filtered 0.01M PBS, pH 7.4 | | |
| Clone: | HENV-26 | | |
| Purity: | > 95% by SDS-PAGE. | | |
| Purification: | Protein A/G affinity chromatography | | |
| Concentration: | 1.0 mg/ml, lot specific | | |
| Applications: | Neutralization, ELISA, FCM | | |
| Storage & Stability: | Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.