

G

Human Anti-HeV/NiV Glycoprotein G (HENV-21) mAb

Catalog No.	RVV08010A RVV08010B	Quantity:	100 µg 1.0 mg
Alternate Names:	Glycoprotein G, Attachment glycoprotein, Receptor binding protein, RBP, Nipah Virus (NiV) / Hendra Virus (HeV)		
Description:	Human neutralizing monoclonal to Hendra virus (isolate Horse/Australia/Hendra/1994) & Nipah virus Glycoprotein G (clone HENV-21)		
UniProt ID (target):	O89343 & Q9IH62		
Specificity:	HeV & NiV Glycoprotein G		
Source:	Mammalian cell line		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	HENV-21		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Neutralization, ELISA, FCM		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com