

## L1

### Human Anti-HPV16 L1/Major capsid protein L1 (H16.U4) mAb

<b>Catalog No.</b>	RVV08804A RVV08804B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Major capsid protein L1, Human Papillomavirus (HPV)		
<b>Description:</b>	Neutralizing human monoclonal to Human papillomavirus type 16 Major capsid protein L1. HPV16, in comparison to HPV types 6 and 11, is more often associated with malignant genital cancers in humans.		
<b>UniProt ID (target):</b>	P03101		
<b>Specificity:</b>	Human papillomavirus type 16 Major capsid protein L1		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG1 kappa		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Clone:</b>	H16.U4		
<b>Purity:</b>	> 95% by SDS-PAGE.		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	Neutralization, ELISA, WB		
<b>Storage &amp; Stability:</b>	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)