

## L1R

### Mouse Anti-Vaccinia virus/VACV Protein L1 (Clone M12B9) mAb

<b>Catalog No.</b>	RVV13902A RVV13902B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Protein L1, Virion membrane protein M25, VACWR088, L1R, Vaccinia virus (VACV)		
<b>Description:</b>	Vaccinia virus is a large, double-stranded DNA virus that carries out its entire replication cycle in the cytoplasm of infected cells. The virus is in the same family as variola virus, the causative agent of smallpox, and immunization with a vaccine derived from vaccinia virus provides protection against smallpox, monkeypox, and other members of the orthopoxvirus family.		
<b>Specificity:</b>	Vaccinia virus Protein L1		
<b>Concentration:</b>	1.0 mg/ml		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Clone:</b>	M12B9		
<b>Purity:</b>	> 95% by SDS-PAGE.		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	Neutralization, ELISA		
<b>Storage &amp; Stability:</b>	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

