

**P/V**

## **Mouse Anti-Parainfluenza virus 5/PIV5 Protein P (Clone SAA0288) mAb**

<b>Catalog No.</b>	RVV14601A RVV14601B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Phosphoprotein, Protein P		
<b>Description:</b>	Neutralizing mouse monoclonal SAA0288 that targets Parainfluenza virus 5 (PIV5) Protein P was purified by affinity chromatography.		
<b>UniProt ID:</b>	P11208		
<b>Specificity:</b>	Parainfluenza virus 5 (PIV5) protein P		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Clone:</b>	SAA0288		
<b>Purity:</b>	> 95% by SDS-PAGE and SEC-HPLC.		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Applications:</b>	Neutralization, ELISA		
<b>Storage &amp; Stability:</b>	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)