

Mouse Anti-Coxsackievirus A6/CVA6 VP1/Capsid protein VP1 (SAA0354) mAb

Catalog No.	RVV19102A RVV19102B	Quantity:	100 µg 1.0 mg
Alternate Names:	Capsid protein VP1, Virion protein 1		
Description:	Mouse Anti-Coxsackievirus A6/CVA6 VP1/Capsid protein VP1 (SAA0354) is an affinity purified, blocking murine monoclonal antibody.		
UniProt ID (target):	A0A0K2BNC7		
Specificity:	Coxsackievirus A6 Capsid protein VP1		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2a		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Functional studies, Blocking, IHC		
Clone:	SAA0354		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

