

gB

Mouse Anti-EBV/HHV4 gB/BALF4 (AMMO5) mAb

Catalog No.	RVV22003A RVV22003B	Quantity:	100 µg 1.0 mg
Alternate Names:	Envelope glycoprotein B, gB, BALF4, gB, EBVaGC1_061, Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)		
Description:	Murine monoclonal to Envelope glycoprotein B from Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4.) The Envelope glycoprotein forms spikes at the surface of virion envelope which is essential for the initial attachment to heparan sulfate moieties of the host cell surface proteoglycans. Involved in fusion of viral and cellular membranes leading to virus entry into the host cell.		
UniProt ID (target):	P0C762		
Specificity:	Envelope glycoprotein B from Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2a kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	AMMO5		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA, FC, Western blot		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com