

Human Anti-HMPV F/Fusion Protein (DS7) mAb

Catalog No.	RVV22902A RVV22902B	Quantity:	100 µg 1.0 mg
Alternate Names:	Fusion glycoprotein F0, F, Protein F		
Description:	Human monoclonal antibody that targets Fusion glycoprotein F0 from human metapneumovirus (HMPV)		
UniProt ID (target):	G3KCK8		
Specificity:	Fusion glycoprotein F0 from human metapneumovirus		
Source:	Mammalian cell line		
Isotype:	Human IgG1 lambda		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	DS7		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA, Neutralization		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com