

C

Anti-HBV-D HBeAg/C/p25 (E6#) mAb

Catalog No.	RVV23601A RVV23601B	Quantity:	100 µg 1.0 mg
Alternate Names:	External core antigen, HBeAg, Precore protein, p25		
Description:	Human monoclonal antibody (clone E6#) specific for the External core antigen from Hepatitis B virus genotype D subtype ayw (isolate France/Tiollais/1979) (HBV-D), expressed in mammalian cells and purified using Protein A or G affinity chromatography.		
UniProt ID (target):	P0C573		
Specificity:	External core antigen from Hepatitis B virus genotype D subtype ayw (isolate France/Tiollais/1979) (HBV-D)		
Source:	Mammalian cell line		
Isotype:	Human IgG1		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	E6#		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com