

GP

Human Anti-MARV GP/Envelope glycoprotein (MR191) mAb

Catalog No.	RVV24301A RVV24301B	Quantity:	100 µg 1.0 mg
Alternate Names:	Envelope glycoprotein, GP1, 2, GP, Virion spike glycoprotein		
Description:	Human monoclonal to Envelope glycoprotein from Lake Victoria marburgvirus (strain Ravn-87) (MARV) (Marburg virus (strainKenya/Ravn/1987)) clone MR191		
UniProt ID (target):	Q1PDC7		
Specificity:	Envelope glycoprotein from Lake Victoria marburgvirus (strain Ravn-87) (MARV) (Marburg virus (strainKenya/Ravn/1987))		
Source:	Mammalian cell line		
Isotype:	Human IgG1		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	MR191		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com