

## Mouse Anti-ASFV Hexon protein p72 (SAb2446) mAb

<b>Catalog No.</b>	RVV48101A RVV48101B RVV48101C	<b>Quantity:</b>	50 µg 100 µg 1.0 mg
<b>Alternate Names:</b>	Hexon protein p72, Major capsid protein, MCP, p72, p73		
<b>Description:</b>	Mouse Anti-ASFV Hexon protein p72 (SAb2446) murine monoclonal antibody that is affinity purified from cell culture supernatant.		
<b>Specificity:</b>	African swine fever virus (ASFV) Hexon protein p72		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cell line		
<b>Isotype:</b>	Mouse IgG		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	SAb2446		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

