


tcdB

Anti-*Clostridioides difficile* Toxin B/tcdB (PA41) mAb

Catalog No.	RXX00303A RXX00303B	Quantity:	100 µg 1.0 mg
Alternate Names:	Toxin B, <i>C. difficile</i>		
Description:	Affinity purified mouse monoclonal antibody that targets Toxin B from <i>Clostridioides difficile</i> (<i>Peptoclostridium difficile</i>)		
UniProt ID:	M4NKT9		
Specificity:	Toxin B from <i>Clostridioides difficile</i> (<i>Peptoclostridium difficile</i>)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2a		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cells		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A or G affinity chromatography from cell culture supernatant		
Clone:	PA41		
Applications:	Neutralization, SPR		
Manufacturer:	Antibody System		
Storage & Stability:	Store at 2-8°C for immediate use, within 2 weeks. Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

SDS-PAGE



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.