

Mouse Anti-LPS/Lipopolysaccharide (SAA0597) mAb

Catalog No.	RXX02126A RXX02126B	Quantity:	100 µg 1.0 mg
Alternate Names:	Lipopolysaccharide (LPS), O-polysaccharide moiety (<i>Klebsiella pneumoniae</i> O3 serotype O-antigen)		
Description:	Recombinant mouse monoclonal antibody which targets Lipopolysaccharide (LPS) from <i>Klebsiella pneumoniae</i> .		
ChEBI ID:	ChEBI: 16412		
Specificity:	Lipopolysaccharide (LPS) from <i>Klebsiella pneumoniae</i> O3 serotype O-antigen		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2a		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cell culture		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Clone:	SAA0597		
Applications:	ELISA, FCM		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

