

fliC

Human Anti-*Escherichia coli* O157:H7 FliC/Flagellin (SAA0364) mAb

Catalog No.	RXX07801A RXX07801B	Quantity:	100 µg 1.0 mg
Alternate Names:	Flagellin, fliC, Z3013, Bacteria, <i>Escherichia coli</i> O157:H7		
Description:	Recombinant human monoclonal antibody which targets Flagellin from <i>Escherichia coli</i> strain O157:H7		
UniProt ID (target):	Q7DBI0		
Source:	Mammalian cells		
Specificity:	Flagellin from <i>Escherichia coli</i> O157:H7		
Isotype:	Human IgG1		
Clone:	SAA0364		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com