

Human Anti-*Staphylococcus aureus* LukAB/LukGH (SAA0520) mAb

Catalog No.	RXX08201A RXX08201B	Quantity:	100 µg 1.0 mg
Alternate Names:	Cytotoxin LukAB, Leukocidin LukGH, LukG, LukA, Leukocidin G, Bi-component leukocidin LukGH subunit G, LukH, LukB, Leukocidin H, Bi-component leukocidin LukGH subunit H, <i>Staphylococcus aureus</i>		
Description:	Recombinant human monoclonal antibody which targets Leukocidin LukAB, LukGH from <i>Staphylococcus aureus</i> strains COL, NCTC 8325 / PS47		
UniProt ID (target):	P21224		
Source:	Mammalian cells		
Specificity:	LukAB/LukGH from <i>Staphylococcus aureus</i> strains COL, NCTC 8325 / PS47		
Isotype:	Human IgG1		
Clone:	SAA0520		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA, FCM		
Storage & Stability:	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.