

Human Anti-Mal f 1/Allergen Pit o 1 (SAA0640) mAb

Catalog No.	RXX08603A RXX08603B	Quantity:	100 µg 1.0 mg
Alternate Names:	Major allergen Mal f 1, Allergen Pit o 1, Mal f 1		
Description:	Recombinant human monoclonal antibody which targets Major allergen Mal f 1 from Malassezia furfur (Pityriasis versicolor infection agent) (Pityrosporum furfur)		
UniProt ID (target):	Q01940		
Source:	Mammalian cells		
Specificity:	Major allergen Mal f 1 from Malassezia furfur (Pityriasis versicolor infection agent) (Pityrosporum furfur)		
Isotype:	VHH-hFc		
Clone:	SAA0640		
Purification:	Protein A or G affinity chromatography from cell culture supernatant		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Concentration:	1.0 mg/ml, lot specific		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

