

fhaB

Human Anti-Bordetella pertussis fhaB/FHA (T5) mAb

Catalog No.	RXX11501A RXX11501B	Quantity:	100 µg 1.0 mg
Alternate Names:	Filamentous hemagglutinin, FHA, fhaB		
Description:	Recombinant human monoclonal antibody which targets filamentous hemagglutinin from <i>Bordetella pertussis</i> (strain Tohama I / ATCC BAA-589 / NCTC 13251)		
UniProt ID:	P12255		
Specificity:	Filamentous hemagglutinin (FHA) from <i>Bordetella pertussis</i> (strain Tohama I / ATCC BAA-589 / NCTC 13251)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Human IgG1		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	XtenCHO cell line		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Clone:	T5		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for 1-2 weeks. Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com