

## Anti-Der p 2 Human IgE (1B8) mAb

<b>Catalog No.</b>	RXX51818A RXX51818B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Description:</b>	Purified human monoclonal IgE that targets Der p 2, Group 2 allergen from Dermatophagoides pteronyssinus (European house dust mite)		
<b>Specificity:</b>	Dermatophagoides pteronyssinus (European house dust mite) Group 2 allergen		
<b>Source:</b>	Mammalian cells		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Isotype:</b>	Human IgE		
<b>Clone:</b>	1B8		
<b>Formulation:</b>	0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Applications:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for 1 week or as aliquots for up to 1 year at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

