

DKK1

InVivoMAb Anti-Human DKK1 (DKN-01)

Catalog No.	VHB48101A VHB48101B	Quantity:	100 µg 1.0 mg
Alternate Names:	Dickkopf-related protein 1, Dickkopf-1, Dkk-1, hDkk-1, SK		
Description:	Mouse monoclonal antibody (clone DKN-01) which targets human Dickkopf-related protein 1		
UniProt ID (target):	O94907		
Source:	Mammalian cells		
Specificity:	Human Dickkopf-related protein 1		
Isotype:	Mouse IgG2a kappa		
Clone:	DKN-01		
Purification:	Protein A or G affinity chromatography from cell culture supernatant		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg by LAL analysis		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Functional studies, Neutralization, ELISA EC ₅₀ = 132.4 ng/ml, binds with human DKK1 by ELISA		
Storage & Stability:	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com