


## IL11


### InVivoMAb Anti-Human IL11 (X203)

<b>Catalog No.</b>	VHD42201A VHD42201B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Interleukin-11, IL-11, Adipogenesis inhibitory factor, AGIF, Oprelvekin		
<b>Description:</b>	Neutralizing murine monoclonal antibody (clone X203) which targets human Interleukin-11		
<b>UniProt ID (target):</b>	P20809		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Human IL-11		
<b>Isotype:</b>	Murine IgG2a kappa		
<b>Clone:</b>	X203		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	Sterile-filtered 100 mM Pro-Ac, 20 mM Arg, pH 5.0		
<b>Endotoxin:</b>	< 0.01 EU/µg		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	Functional studies, Neutralizing		
<b>Application Notes:</b>	EC <sub>50</sub> = 16.39 ng/ml, clone X203 binds with human IL-11 by indirect ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

SEC-HPLC



Anti-Human IL11 Antibody (X203) binds with IL11 protein



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)