

IL11RA

InVivoMAb Anti-Human IL11RA (X209)

Catalog No.	VHG84901A VHG84901B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-11 receptor subunit alpha, IL-11 receptor subunit alpha, IL-11R subunit alpha, IL-11R-alpha, IL-11RA		
Description:	Neutralizing human monoclonal antibody (clone X209) which targets Interleukin-11 receptor subunit alpha.		
UniProt ID (target):	Q14626		
Source:	Mammalian cells		
Specificity:	Human Interleukin-11 receptor subunit alpha		
Isotype:	Mouse IgG2a kappa		
Clone:	X209		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Functional studies, Neutralizing, ELISA		
Storage & Stability:	Store at 2-8°C for immediate use/1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com