

## Cd19 & Cd3e

### InVivoMAb Anti-Mouse CD19 & CD3e Bispecific Antibody (Iv0224)

<b>Catalog No.</b>	VMD10801A VMD10801B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	B-lymphocyte antigen CD19, Differentiation antigen CD19, CD19, T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e		
<b>Description:</b>	Bispecific murine monoclonal antibody (clone Iv0224) which targets mouse CD19 and CD3e purified with low endotoxin.		
<b>UniProt ID (target):</b>	P25918 & P22646		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Mouse CD19 and CD3e		
<b>Isotype:</b>	Murine IgG2a kappa		
<b>Clone:</b>	Iv0224		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	Sterile-filtered PBS, pH 6.0 containing 0.2M Arg <b>Note: formulation change</b>		
<b>Endotoxin:</b>	< 0.01 EU/µg by LAL analysis		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)