

IL17a

InVivoMAb Anti-Mouse IL17A (BZN035)

Catalog No.	VMH28801A VMH28801B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-17A, IL-17A, Cytotoxic T-lymphocyte-associated antigen 8, CTLA-8		
Description:	Neutralizing mouse monoclonal antibody (clone BZN035) which targets mouse Interleukin-17A		
UniProt ID (target):	Q62386		
Source:	Mammalian cells		
Specificity:	Mouse IL-17A		
Isotype:	Mouse IgG2a		
Clone:	BZN035		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Functional studies, Neutralizing, FCM		
Storage & Stability:	Store at 2-8°C for immediate use/1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.