

## Trem2

### InVivoMAb Anti-Mouse TREM2 (4D9)

<b>Catalog No.</b>	VMJ69201A VMJ69201B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Triggering receptor expressed on myeloid cells 2, TREM-2, Triggering receptor expressed on monocytes 2		
<b>Description:</b>	Affinity purified, low endotoxin murine monoclonal anti-Mouse TREM2 may be used as an agonist in functional studies.		
<b>UniProt ID (target):</b>	Q99NH8		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	Mouse TREM-2		
<b>Isotype:</b>	Mouse IgG2a		
<b>Clone:</b>	4D9		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Endotoxin:</b>	< 1 EU/mg		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Applications:</b>	Functional studies (agonist), ELISA, FCM		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use within 1-2 weeks or aliquot the stock solution and store at -20°C to -80°C up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

